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Designed for applying large-intestinal and large-and-small intestinal circular compression anastomoses on the rectum and colon which differs essentially from the available staple devices.

Anastomoses are applied by compressing the tissue of the organs to be connected, with a subsequent natural separation of the compressing elements 7 to 14 days after the operation.

The edges of the hollow organs on which an anastomosis is applied are connected throughout their circumference by compressing elements when the handle of the device is pressed.

Simultaneously with a circular suture, the device dissects a surplus tissue between the organs sutured, the dissected tissue of the intestinal walls being removed together with the withdrawal of the device.

The device is supplied complete with an appliance which improves the aseptic safety of the operation and obviates the need for application of a purse-string suture on the distal part of the intestine before the anastomosis.

The АКА-2 device operates in the same way as the staple devices ПКС, СК-28, СПТУ, therefore manipulation of the device does not require any special skill and can be used by any surgeon.

The radical difference of the АКА-2 from other staple devices provides for the following advantages:

- after necrotization of intestinal wall portions compressed by the connecting rings, the latter are detached, move into the large intestine lumen and excreted naturally without leaving any foreign bodies in the anastomosis area or causing any complications;

- the anastomosis is highly tight and reliable, its healing proceeds through primary tensioning of the tissues;


- the operation time is reduced, because the second row of manual suture is not required due to the proper adjustment of the intestine ends to be sutured, and the absence of tension;

- the device eliminates the possibility of constriction or formation of a scar constriction in the anastomosis area;

- better aseptic conditions are provided which minimizes the danger of wound suppuration and prevents failure of the anastomosis;

- there is no need in applying a purse-string suture which is particularly difficult when the rectum stump has a low position;

- it is much easier to check the condition of the device before the operation; the time of preparation of the device for operation and the operation time

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